

# SOP pre-processing MethyLight raw data qWID-EC

Last updated: 10.01.2023, Charlotte Vavourakis

1. **Save the .eds file in the ZIDshare**
2. **Make a copy of the .eds file to your project folder and experiment subfolder in:**

\Dropbox\eutops\data\raw\_data\methyLight\

3. **Generate a samplesheet using the Shiny app**

- Prepare a sample sheet in Excel, first column should be labeled "Number", second column should be labeled "Sample\_name". Up to 90 samples can be entered, mind the trailing "0" for the first 9 samples! Example:

	A	B	C
1	Number	Sample_name	
2	Sample_01	test1	
3	Sample_02	test2	
4	Sample_03	test3	
5	Sample_04	test4	
6	Sample_05	test5	
7	Sample_06	test6	
8	Sample_07	test7	
9	Sample_08	test8	
10	Sample_09	test9	
11	Sample_10	test10	
12	Sample_11	test11	
13	Sample_12	test12	
14	Sample_13	test13	

- Make sure R is installed and in your environment path. Libraries needed: shiny, readxl, stringr and dplyr.
- Navigate to \Dropbox\eutops\scripts\methyLight\ML\_WIDqGyn\shiny
- For Windows: Double-click GeneratePlatelayouEC.bat (you could make a shortcut to this file e.g. on your Desktop)
- Follow instructions in the browser.
- Please note that .bat files are specific for Windows OS.
- Note, your samples should be pipetted like this


	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	EA-0001 25CARN2 GVPC1	EA-0019 25CARN2 GVPC1	EA-0027 25CARN2 GVPC1	EA-0028 25CARN2 GVPC1	EA-0043 25CARN2 GVPC1	EA-0051 25CARN2 GVPC1	EA-0059 25CARN2 GVPC1	EA-0067 25CARN2 GVPC1	EA-0075 25CARN2 GVPC1	EA-0085 25CARN2 GVPC1	Sample X 25CARN2 GVPC2	EA-0001 COL2A1	EA-0009 COL2A1	EA-0019 COL2A1	EA-0027 COL2A1	EA-0028 COL2A1	EA-0043 COL2A1	EA-0051 COL2A1	EA-0059 COL2A1	EA-0067 COL2A1	EA-0075 COL2A1	EA-0085 COL2A1	Sample X COL2A1	
B	EA-0001 25CARN2 GVPC1	EA-0009 25CARN2 GVPC1	EA-0019 25CARN2 GVPC1	EA-0027 25CARN2 GVPC1	EA-0028 25CARN2 GVPC1	EA-0043 25CARN2 GVPC1	EA-0051 25CARN2 GVPC1	EA-0059 25CARN2 GVPC1	EA-0067 25CARN2 GVPC1	EA-0075 25CARN2 GVPC1	Sample X 25CARN2 GVPC2	EA-0001 COL2A1	EA-0009 COL2A1	EA-0019 COL2A1	EA-0027 COL2A1	EA-0028 COL2A1	EA-0043 COL2A1	EA-0051 COL2A1	EA-0059 COL2A1	EA-0067 COL2A1	EA-0075 COL2A1	EA-0085 COL2A1	Sample X COL2A1	
C	EA-0002 25CARN2 GVPC1	EA-0001 25CARN2 GVPC1	EA-0029 25CARN2 GVPC1	EA-0028 25CARN2 GVPC1	EA-0044 25CARN2 GVPC1	EA-0052 25CARN2 GVPC1	EA-0060 25CARN2 GVPC1	EA-0069 25CARN2 GVPC1	EA-0077 25CARN2 GVPC1	EA-0086 25CARN2 GVPC1	Sample Y 25CARN2 GVPC2	EA-0002 COL2A1	EA-0001 COL2A1	EA-0029 COL2A1	EA-0028 COL2A1	EA-0044 COL2A1	EA-0052 COL2A1	EA-0060 COL2A1	EA-0069 COL2A1	EA-0077 COL2A1	EA-0086 COL2A1	Sample Y COL2A1		
D	EA-0002 25CARN2 GVPC1	EA-0001 25CARN2 GVPC1	EA-0029 25CARN2 GVPC1	EA-0028 25CARN2 GVPC1	EA-0044 25CARN2 GVPC1	EA-0052 25CARN2 GVPC1	EA-0060 25CARN2 GVPC1	EA-0069 25CARN2 GVPC1	EA-0077 25CARN2 GVPC1	EA-0086 25CARN2 GVPC1	Sample Y 25CARN2 GVPC2	EA-0002 COL2A1	EA-0001 COL2A1	EA-0029 COL2A1	EA-0028 COL2A1	EA-0044 COL2A1	EA-0052 COL2A1	EA-0060 COL2A1	EA-0069 COL2A1	EA-0077 COL2A1	EA-0086 COL2A1	Sample Y COL2A1		
E	EA-0003 25CARN2 GVPC1	EA-0002 25CARN2 GVPC1	EA-0021 25CARN2 GVPC1	EA-0029 25CARN2 GVPC1	EA-0027 25CARN2 GVPC1	EA-0045 25CARN2 GVPC1	EA-0053 25CARN2 GVPC1	EA-0061 25CARN2 GVPC1	EA-0069 25CARN2 GVPC1	EA-0078 25CARN2 GVPC1	EA-0087 25CARN2 GVPC1	STD1 25CARN2 GVPC2	EA-0003 COL2A1	EA-0002 COL2A1	EA-0021 COL2A1	EA-0029 COL2A1	EA-0027 COL2A1	EA-0045 COL2A1	EA-0053 COL2A1	EA-0061 COL2A1	EA-0069 COL2A1	EA-0078 COL2A1	EA-0087 COL2A1	STD1 COL2A1
F	EA-0003 25CARN2 GVPC1	EA-0002 25CARN2 GVPC1	EA-0021 25CARN2 GVPC1	EA-0029 25CARN2 GVPC1	EA-0027 25CARN2 GVPC1	EA-0045 25CARN2 GVPC1	EA-0053 25CARN2 GVPC1	EA-0061 25CARN2 GVPC1	EA-0069 25CARN2 GVPC1	EA-0078 25CARN2 GVPC1	EA-0087 25CARN2 GVPC1	STD1 25CARN2 GVPC2	EA-0003 COL2A1	EA-0002 COL2A1	EA-0021 COL2A1	EA-0029 COL2A1	EA-0027 COL2A1	EA-0045 COL2A1	EA-0053 COL2A1	EA-0061 COL2A1	EA-0069 COL2A1	EA-0078 COL2A1	EA-0087 COL2A1	STD1 COL2A1
G	EA-0004 25CARN2 GVPC1	EA-0003 25CARN2 GVPC1	EA-0022 25CARN2 GVPC1	EA-0030 25CARN2 GVPC1	EA-0038 25CARN2 GVPC1	EA-0046 25CARN2 GVPC1	EA-0054 25CARN2 GVPC1	EA-0062 25CARN2 GVPC1	EA-0070 25CARN2 GVPC1	EA-0079 25CARN2 GVPC1	STD2 25CARN2 GVPC2	EA-0004 COL2A1	EA-0003 COL2A1	EA-0022 COL2A1	EA-0030 COL2A1	EA-0038 COL2A1	EA-0046 COL2A1	EA-0054 COL2A1	EA-0062 COL2A1	EA-0070 COL2A1	EA-0079 COL2A1	EA-0088 COL2A1	STD2 COL2A1	
H	EA-0004 25CARN2 GVPC1	EA-0003 25CARN2 GVPC1	EA-0022 25CARN2 GVPC1	EA-0030 25CARN2 GVPC1	EA-0038 25CARN2 GVPC1	EA-0046 25CARN2 GVPC1	EA-0054 25CARN2 GVPC1	EA-0062 25CARN2 GVPC1	EA-0070 25CARN2 GVPC1	EA-0079 25CARN2 GVPC1	STD2 25CARN2 GVPC2	EA-0004 COL2A1	EA-0003 COL2A1	EA-0022 COL2A1	EA-0030 COL2A1	EA-0038 COL2A1	EA-0046 COL2A1	EA-0054 COL2A1	EA-0062 COL2A1	EA-0070 COL2A1	EA-0079 COL2A1	EA-0088 COL2A1	STD2 COL2A1	
I	EA-0005 25CARN2 GVPC1	EA-0004 25CARN2 GVPC1	EA-0023 25CARN2 GVPC1	EA-0031 25CARN2 GVPC1	EA-0039 25CARN2 GVPC1	EA-0047 25CARN2 GVPC1	EA-0055 25CARN2 GVPC1	EA-0063 25CARN2 GVPC1	EA-0071 25CARN2 GVPC1	EA-0080 25CARN2 GVPC1	STD3 25CARN2 GVPC2	EA-0005 COL2A1	EA-0004 COL2A1	EA-0023 COL2A1	EA-0031 COL2A1	EA-0039 COL2A1	EA-0047 COL2A1	EA-0055 COL2A1	EA-0063 COL2A1	EA-0071 COL2A1	EA-0080 COL2A1	EA-0089 COL2A1	STD3 COL2A1	
J	EA-0005 25CARN2 GVPC1	EA-0004 25CARN2 GVPC1	EA-0023 25CARN2 GVPC1	EA-0031 25CARN2 GVPC1	EA-0039 25CARN2 GVPC1	EA-0047 25CARN2 GVPC1	EA-0055 25CARN2 GVPC1	EA-0063 25CARN2 GVPC1	EA-0071 25CARN2 GVPC1	EA-0080 25CARN2 GVPC1	STD3 25CARN2 GVPC2	EA-0005 COL2A1	EA-0004 COL2A1	EA-0023 COL2A1	EA-0031 COL2A1	EA-0039 COL2A1	EA-0047 COL2A1	EA-0055 COL2A1	EA-0063 COL2A1	EA-0071 COL2A1	EA-0080 COL2A1	EA-0089 COL2A1	STD3 COL2A1	
K	EA-0006 25CARN2 GVPC1	EA-0005 25CARN2 GVPC1	EA-0024 25CARN2 GVPC1	EA-0032 25CARN2 GVPC1	EA-0040 25CARN2 GVPC1	EA-0048 25CARN2 GVPC1	EA-0056 25CARN2 GVPC1	EA-0064 25CARN2 GVPC1	EA-0072 25CARN2 GVPC1	EA-0081 25CARN2 GVPC1	STD4 25CARN2 GVPC2	EA-0006 COL2A1	EA-0005 COL2A1	EA-0024 COL2A1	EA-0032 COL2A1	EA-0040 COL2A1	EA-0048 COL2A1	EA-0056 COL2A1	EA-0064 COL2A1	EA-0072 COL2A1	EA-0081 COL2A1	EA-0090 COL2A1	STD4 COL2A1	
L	EA-0006 25CARN2 GVPC1	EA-0005 25CARN2 GVPC1	EA-0024 25CARN2 GVPC1	EA-0032 25CARN2 GVPC1	EA-0040 25CARN2 GVPC1	EA-0048 25CARN2 GVPC1	EA-0056 25CARN2 GVPC1	EA-0064 25CARN2 GVPC1	EA-0072 25CARN2 GVPC1	EA-0081 25CARN2 GVPC1	STD4 25CARN2 GVPC2	EA-0006 COL2A1	EA-0005 COL2A1	EA-0024 COL2A1	EA-0032 COL2A1	EA-0040 COL2A1	EA-0048 COL2A1	EA-0056 COL2A1	EA-0064 COL2A1	EA-0072 COL2A1	EA-0081 COL2A1	EA-0090 COL2A1	STD4 COL2A1	
M	EA-0007 25CARN2 GVPC1	EA-0006 25CARN2 GVPC1	EA-0025 25CARN2 GVPC1	EA-0033 25CARN2 GVPC1	EA-0041 25CARN2 GVPC1	EA-0049 25CARN2 GVPC1	EA-0057 25CARN2 GVPC1	EA-0065 25CARN2 GVPC1	EA-0073 25CARN2 GVPC1	EA-0082 25CARN2 GVPC1	STD5 25CARN2 GVPC2	EA-0007 COL2A1	EA-0006 COL2A1	EA-0025 COL2A1	EA-0033 COL2A1	EA-0041 COL2A1	EA-0049 COL2A1	EA-0057 COL2A1	EA-0065 COL2A1	EA-0073 COL2A1	EA-0082 COL2A1	EA-0091 COL2A1	STD5 COL2A1	
N	EA-0007 25CARN2 GVPC1	EA-0006 25CARN2 GVPC1	EA-0025 25CARN2 GVPC1	EA-0033 25CARN2 GVPC1	EA-0041 25CARN2 GVPC1	EA-0049 25CARN2 GVPC1	EA-0057 25CARN2 GVPC1	EA-0065 25CARN2 GVPC1	EA-0073 25CARN2 GVPC1	EA-0082 25CARN2 GVPC1	STD5 25CARN2 GVPC2	EA-0007 COL2A1	EA-0006 COL2A1	EA-0025 COL2A1	EA-0033 COL2A1	EA-0041 COL2A1	EA-0049 COL2A1	EA-0057 COL2A1	EA-0065 COL2A1	EA-0073 COL2A1	EA-0082 COL2A1	EA-0091 COL2A1	STD5 COL2A1	
O	EA-0008 25CARN2 GVPC1	EA-0007 25CARN2 GVPC1	EA-0026 25CARN2 GVPC1	EA-0034 25CARN2 GVPC1	EA-0042 25CARN2 GVPC1	EA-0050 25CARN2 GVPC1	EA-0058 25CARN2 GVPC1	EA-0066 25CARN2 GVPC1	EA-0074 25CARN2 GVPC1	EA-0084 25CARN2 GVPC1	H2O 25CARN2 GVPC2	EA-0008 COL2A1	EA-0007 COL2A1	EA-0026 COL2A1	EA-0034 COL2A1	EA-0042 COL2A1	EA-0050 COL2A1	EA-0058 COL2A1	EA-0066 COL2A1	EA-0074 COL2A1	EA-0084 COL2A1	EA-0092 COL2A1	H2O COL2A1	
P	EA-0008 25CARN2 GVPC1	EA-0007 25CARN2 GVPC1	EA-0026 25CARN2 GVPC1	EA-0034 25CARN2 GVPC1	EA-0042 25CARN2 GVPC1	EA-0050 25CARN2 GVPC1	EA-0058 25CARN2 GVPC1	EA-0066 25CARN2 GVPC1	EA-0074 25CARN2 GVPC1	EA-0084 25CARN2 GVPC1	H2O 25CARN2 GVPC2	EA-0008 COL2A1	EA-0007 COL2A1	EA-0026 COL2A1	EA-0034 COL2A1	EA-0042 COL2A1	EA-0050 COL2A1	EA-0058 COL2A1	EA-0066 COL2A1	EA-0074 COL2A1	EA-0084 COL2A1	EA-0092 COL2A1	H2O COL2A1	

#### 4. Save the resulting sample sheet, alongside the raw data (.eds file) both on the ZIDshare and the Dropbox


!!!!!!! Make sure to do this BEFORE editing an .eds file !!!!!!!

#### 5. Using the Quantstudio Design & Analysis software 2.6.0, extract, normalize and export run results

- Open .eds file
- Import your sample sheet under Plate setup > ... > Import Plate Setup



- Press ok, when prompted to overwrite Plate Setup



- Under Plate setup, Add the correct Quantities for the standards in the Sample list (example given for gBlocks)

Design & Analysis Software 2.0.0

DA Data Gallery 2022-01-27\_apoptotic-check\_20220127\_135923.edb

Run Method Plate Setup Run Summary Quality Check Standard Curve

Target SNP Passive Reference ROX View

Well	Sample Name	Sample Type	Target Name	Target Reporter	Target Quencher
A1					
A2					
A3					
A4					
A5					
A6					
A7					
A8					
A9					
A10					
A11					
A12					
A13					
A14					
A15					
A16					
A17					
A18					

**Samples (16)** Biogroup (3)


Name	Color	Type	Quantity
Nalm6_apo_E_BM-D...	Yellow	Unknown	
Nalm6_apo_L_BM-D...	Yellow	Unknown	
Nalm6_ctl_BM-DNA	Yellow	Unknown	
STD1	Blue	Standard	1,000,000
STD2	Red	Standard	100,000
STD3	Green	Standard	10,000
STD4	Red	Standard	1,000
gBlock	Purple	Unknown	

**Targets (3)** Reagents (3)


Target	Color	Reporter	Quencher
COL3A1	Red	FAM	NFQ-MGB
EpC	Green	CY5	NFQ-MGB
hMC	Yellow	VIC	NFQ-MGB

- Press Analyze
- Under Quality check inspect the Amplification Plots
- You can select specific wells, checking minimally Standards, Positive and Negative controls.

- Delta Rn should be normalized on the Y-axis using the passive Reference dye ROX, showing a scale from 0 to 1



- Under Actions, Primary Analysis Settings, Set the Thresholds, and press Save:
  - COL2A1 (FAM) = 0,12
  - ZSCAN12 (FAM) = 0.17
  - GYPC1 (CY5), GYPC2 (CY5) = 0.04



- Under Actions export the result files (default setting)
- Save them to your project folder and experiment subfolder in:

\\Dropbox\\eutops\\data\\raw\_data\\methyLight\\

- Close and save the .eds file

**6. Calculate PMR, WID-qEC and WID-qCIN using Shiny**

- Navigate to \Dropbox\eutops\scripts\methyLight\ML\_WIDqGyn\shiny
- For Windows: double-click CalculatePMRGyn.bat (you could make a shortcut to this file e.g. on your Desktop)
- For Mac/Linux: double-click CalculatePMRGyn.sh (you could make a shortcut to this file e.g. on your Desktop)
- Follow instructions in the browser.

**7. Please update “experiment\_log” \Dropbox\eutops\data! For larger cohorts and studies**